## **MOUNTING PRESSES** (ADVANCED TYPE)





- Built-in heating and cooling structure reduces the mounting time to 6 minutes
  The programmed control system reserves 10 customized preparation schemes to meet different mounting requirements
- The thickness of the middle mold is up to 15mm, and 2 samples can be mounted at the same time at one station
- Can accurately control the pressure size, adjustment range: 1~20MPa



hot mounting resin(included)

## **SPECIFICATION**

Code	MLP-MP43	MLP-MP53	
Station number	1	2	
Power	600W+Heating power*	700W+Heating power×2	
Mounting diameter	Ø30mm(Customizable Ø25, Ø40, Ø50mm)		
System pressure	1~20MPa, Optional increase or decrease 1MPa		
Pressure mode	Automatic electro-hydraulic		
Cooling mode	Automatic water cooling**		
Heating range	90~200°C		
Holding time	0~10min		
Cooling temperature/time	40~90°C/0~10min		
Power supply	AC 220V, 50Hz		
Weight	44kg	77kg	
Dimension	530×485×463mm	671×485×463mm	

<sup>\*</sup>Heating power of mounting sleeve: 1800W (mounting diameter≤31.8mm), 2400W(mounting diameter≥38.1mm)

## STANDARD DELIVERY

STANDARD DELIVERT		
Main unit	1рс	
Hot mounting resin	1 bottle (250g)	
Funnel	1pc	
Filter	1pc	
vacuum cleaner tip	1pc	

## **OPTIONAL ACCESSORY**

Description	Code	Specification
Hot mounting resin	MLP-MP30-PDR	black,4kg

<sup>\*\*\*</sup>Cooling water pressure should not be higher than 1MPa, must be filtered through 40µm filter